

Things to consider when selecting a scooter or power chair:

Who is using the power mobility device?

- What does the person weigh and how tall is he or she?
- Does the person have any physical limitations?

Where are you going to be using the scooter or power chair?

- Indoors, outdoors or both?
- If outdoors, on what type of surface will you be riding?

Do you want a 3 wheel or 4 wheel model?

- 3 Wheel Scooters have a tighter turning radius than a 4 Wheel (good for tighter spaces). Once broken down for transport, a 3 Wheel Scooter is lighter to pick up than a 4 Wheel Scooter and they have more leg room. 4 Wheel Scooters are great for outdoor riding especially around in a yard and over rougher terrain.

If the mobility device will be used away from your home, how will you transport it?

- Do you want to transport it inside your vehicle?
 - Are you looking for a portable model that comes apart easily?
 - If so how much can you lift?
 - Is a vehicle lift or a portable ramp an better choice for you?
- Is an outside mounted vehicle lift the best choice for you?

Avenue Medical's friendly and experienced staff will explain the different options and ask all the right questions to make sure that you select the best scooter or power wheelchair for you. There are many things to consider but we'll make it easy!