



COMMODES

Model #: _____ Serial #: _____

Setting up Your Equipment:

Make sure commode is adjusted to proper height for user.

Using Your Equipment:

- Seat must be down before using commode; do not sit on pail holder. (*Diagram 1*)
- If using a mobile commode, be sure wheels are locked before using.
- If commode has drop-arms, ensure they are locked in the upright position before using. (*Diagram 2*)
- If using elevated commode seat on toilet bowl, make sure the seat is properly attached to prevent loss of balance and/or tipping.

Diagram 1




Diagram 2



Maintenance:

- Make sure all screws; nuts and bolts are tight at all times.
- Inspect rubber tips on the legs to ensure they are not worn or torn.
- For hygiene sake, empty pail and clean it after each use.
- Make sure wheels and locks are working properly.

Safety Issues:

- Do not sit on commode for extended periods of time.
- Commode seats are not intended to serve as wheelchairs.
- Make sure the user can support themselves' while using the commode. If user is unsteady or weak, a caregiver should provide support during the toileting process.

Available Accessories (*not all accessories are listed here – please ask staff about other items*):

*Cleaning
Wipes*



*Replacement
Pails*



*Splash
Guards*



*Disposable
Liner Bags*

